

Recover more volume. Recover it faster. Filter sterilize with the best.

Introducing the NEW and improved Millipore Express® PLUS membrane for Stericup® and Steriflip® filters to ensure fast and complete recovery, with no loss of media or precious additives in the filtration process. With the best sterile filters, you'll do your best work.

SPECIAL OFFER

Buy 3 cases in any combination of Stericup®, Steritop® & Steriflip® filters

GET ONE CASE FREE!*

See reverse for offer details & ordering information
To learn more about sterile filtration products, visit: www.fishersci.com/emdmillipore


Offer expires December 31, 2014.




Still the best.

Stericup® filters with the all NEW Millipore Express® PLUS PES membrane

Faster filtration of media (DMEM/10% FBS, 500 mL)



Least nonspecific binding (radiolabeled IgG)



Less hold up volume: With a membrane structure that doesn't trap your filtrate, you'll recover more of your precious media and supplements when you sterilize them with Stericup® filters.

Stericup®, Steritop® and Steriflip® Filters: The Industry Standard for Sterile Filtration Quality and Efficiency

Stericup® Filter Units – combine a filter unit with a receiver flask and cap for processing and storage.

Description	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Receiver Bottle (mL)	Qty/Pk	Catalogue No.
Stericup®-VP Filter Units	Millipore Express® (PES)	0.1	250	250	12	SCVPU02RE
			1000	1000	12	SCVPU11RE
Our featured 250 mL and 1000 mL Stericup® units provide fast flow and low-protein-binding 0.1 µm PES membrane for sterilization and mycoplasma removal!						
Stericup®-GP Filter Units STEM CELL TESTED	Millipore Express® PLUS (PES)	0.22	150	150	12	SCGPU01RE
			250	250	12	SCGPU02RE
			500	500	12	SCGPU05RE
			500	1000	12	SCGPU10RE
			1000	1000	12	SCGPU11RE
Stericup®-GV Filter Units	Durapore® (PVDF)	0.22	150	150	12	SCGVU01RE
			250	250	12	SCGVU02RE
			500	500	12	SCGVU05RE
			500	1000	12	SCGVU10RE
			1000	1000	12	SCGVU11RE
Stericup®-HV Filter Units	Durapore® (PVDF)	0.45	150	150	12	SCHVU01RE
			250	250	12	SCHVU02RE
			500	500	12	SCHVU05RE
			1000	1000	12	SCHVU11RE

Steritop® Bottle Top Filter Units – can be used with bottles having a 33 mm or 45 mm opening.

Description	Membrane	Pore Size (µm)	Funnel Capacity (mL)	Receiver Bottle (mL)	Qty/Pk	Catalogue No.
Steritop®-GP Filter Units STEM CELL TESTED	Millipore Express® PLUS (PES)	0.22	150	33	12	SCGPS01RE
			150	45	12	SCGPT01RE
			250	33	12	SCGPS02RE
			250	45	12	SCGPT02RE
			500	33	12	SCGPS05RE
			500	45	12	SCGPT05RE
			1000	45	12	SCGPT10RE
Steritop®-GV Filter Units	Durapore® (PVDF)	0.22	500	45	12	SCGVT05RE

Steriflip® Filter Units – completely eliminate liquid handling for sterile samples you can trust.

Description	Catalogue No.
0.22 µm Millipore Express® (PES) Steriflip®-GP Filter unit	SCGP00525
0.45 µm Durapore® (PVDF) Steriflip®-HV Filter unit	SE1M003M00

To receive your free case of filters, simply fax the completed form along with your proof of purchase* to 1-781-533-3280 or E-mail to OSM@merckgroup.com by December 31, 2014. Please allow 6 weeks for processing.

NAME _____

COMPANY _____

ADDRESS _____

CITY _____

STATE/PROVIDENCE _____ ZIP/POSTAL CODE _____

PHONE _____ EMAIL _____

CATALOGUE NO. OF REQUESTED FREE PRODUCT _____

Proof of purchase consists of one of the following: Order confirmation, delivery ticket, packing list, order summary, invoice or acknowledgment. PA121913TM
Limit: one order per customer. This promotion is for one-time-use only and cannot be duplicated or transferred.



To Receive Technical Assistance
in the U.S. and Canada

call toll free: 1(800)645-5476

For Technical Service, please email to:
www.emdmillipore.com/techservice



**Fisher
Scientific**

Part of Thermo Fisher Scientific

In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: www.fishersci.com

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: www.fishersci.ca

www.fishersci.com/emdmillipore