



Multichannel CappAero96 multichannel Pipette & Channels

Easily exchangeable counter

• Choose between variable or fixed volume. Lockable variable volume controllers. Fixed controllers for routine work available on request.

Easy tip ejection

• Large pad for easy tip ejection. Colorcoded acc. to pipette's volume range.

Fully autoclavable body

• Recommended procedure: 20 min. at 121°C -no disassembly required.

Robust design & ultimate precision

• Metal tip cones with double O-rings for optimal sealing and precision. Compatible with all major pipette tip brands. Low incidence of breakage resulting from accidental dropping or tips mounting.

Cat# C10-8, C100-8, C200-8, and C300-8



Contact Van for details:

(647) 919-5154

Van@froggabio.com

