


The ONE You Have Been Waiting For



OneTaq™ DNA Polymerase OneTaq™ Hot Start DNA Polymerase

Why use more than one polymerase for routine and difficult PCR experiments? With OneTaq DNA Polymerase, achieve robust yields with minimal optimization across a wide range of templates. For routine, AT- or GC-rich amplicons, make OneTaq DNA Polymerase your first choice for endpoint PCR.

- Robust yield with minimal optimization
- Ideal for routine, AT- or GC-rich templates
- Hot start, master mix and Quick Load versions available
- Compatible with standard Taq protocols



Amplification of a selection of sequences with varying AT and GC content from human and *C. elegans* genomic DNA, using OneTaq DNA Polymerase.

| | | MedStore Price |
|--------|---|----------------|
| M0480L | OneTaq™ DNA Polymerase - 1000 units | \$186.56 |
| M0480S | OneTaq™ DNA Polymerase - 200 units | \$46.64 |
| M0480X | OneTaq™ DNA Polymerase - 5000 units | \$748.88 |
| M0481L | OneTaq™ Hot Start DNA Polymerase - 1000 units | \$388.96 |
| M0481S | OneTaq™ Hot Start DNA Polymerase - 200 units | \$96.80 |
| M0481X | OneTaq™ Hot Start DNA Polymerase - 5000 units | \$1,556.72 |

New England Biolabs, Ltd. Telephone:
1-800-387-1095; Fax: 1-800-563-3789; email: info@ca.neb.com
Exclusive Canadian distributor for Cell Signaling Technology, Inc.

All online orders receive a
10% discount and free shipping

For more information, please visit www.neb.ca

