

GoTaq[®] G2 Hot Start Master Mixes

Enjoy the Convenience of Room Temperature, Master Mix PCR Setup!

Ready-to-use mixes that contain all necessary components, including the 2nd-generation GoTaq G2 Hot Start Polymerase, for your PCR. Add template and go.

- Simplify reaction setup and save time with a ready-to-use master mix
- Prepare your reaction at room temperature, not on ice
- Eliminate nonspecific amplification with hot start enzyme
- Choose from green or colorless options
- Use at no risk: Backed by the Promega PCR Satisfaction Guarantee

Thursday, December 3rd

10:00am- 12:00pm

University of Toronto MedStore

Free Coffee & Food, too!

Ask Lucy for details...

Lucy Lu

Client Support Associate

lucy.lu@promega.com

514.518.0496

**Buy Any Promega Product at
the Show and Get a**

**FREE GoTaq[®] G2 Hot Start
Green Master Mix Sample!**

GoTaq is a registered trademark of Promega Corporation. Products may be covered by issued or pending patents or subject to use limitations, please visit www.promega.com for more information.

