

Ready-to-use, one step protein gel staining solutions!

- *Fast (5-60 min) protein staining in a single step*
- *No fixing or washing required*
- *No hazardous methanol or acetic acid, stains are entirely aqueous-based*
- *Fully compatible with mass spectrometry and Edman-based sequencing*

One-Step Blue™ stained proteins are detected by visible blue staining or by near-infrared fluorescence.

Cat Number: 21003-1L
Price: \$125.72 CAD per litre

One-Step Lumitein™ is a sensitive red fluorescent gel stain, for protein detection using a UV transilluminator or laser gel scanner.

Cat Number: 21004-1L
Price: \$161.00 CAD per litre