



## Invitrogen Pre-Cast Protein Gel Promotion through MedStore!

For a full list of our pre-cast protein gels please visit

[www.thermofisher.com/proteingels](http://www.thermofisher.com/proteingels)

Promo valid until June 29<sup>th</sup>, 2018.

**NOVEX** WedgeWell Tris-Glycine Gels - **\$79** per box of 10 gels

**BOLT** Bis-Tris Plus Gels - **\$89** per box of 10 gels

**NUPAGE** Bis-Tris Gels- **\$149** per box of 10 gels

**NUPAGE** Tris-Acetate Gels- **\$149** per box of 10 gels



For any questions or concerns, please contact:

Sanam Taheri | Biosciences Account Manager | [Sanam.Taheri@thermofisher.com](mailto:Sanam.Taheri@thermofisher.com)

Patrick Policelli | Supply Center Representative | [Patrick.Policelli@thermofisher.com](mailto:Patrick.Policelli@thermofisher.com)

## Interested in streamlining your Western blotting process? Check out these instruments:

### Power Blotter XL Transfer Device



- Transferring low MW, mixed range MW and high MW samples in 5-10 minutes.
- Run up to 4 mini-sized gels or 2 midi-sized gels at one time.
- Option to use more convenient Power Blotter Select Stacks for pre-assembled transfer stacks.

### iBind Flex Detection System



- The iBind Flex is an automated detection system using sequential lateral flow with a run time of less than 3 hours.
- The iBind Flex is compatible with all of your chromogenic, chemiluminescent, or fluorescent western blotting protocols.
- Run up to one midi, two mini, or 6 vertical strips.

### iBright Imaging System



- The iBright Imaging system captures images and analyzes data from western blots, protein gels and nucleic acid gels.
- Includes a 9.1MP high resolution camera.
- Ability to achieve a 4-plex fluorescent western.
- Advanced functionality with auto sample alignment, autofocus, zoom, and exposure.
- Can image up to 4 blots simultaneously for higher throughput.

Please contact your Biosciences Account Manager if you're interested in a demo:

Sanam Taheri | Biosciences Account Manager | [Sanam.Taheri@thermofisher.com](mailto:Sanam.Taheri@thermofisher.com)