

DNA Ladder **SALE!**


BUY ANY 10

**100bp or 1kb DNA Ladder Plus
50µg & RECEIVE**

5 FREE

100bp ladder plus

M1071/M1072




5 µl/lane, 8 cm gel
0.5×TBE, 5 V/cm, 1h

Advised Electrophoresis Conditions

1kb ladder plus

M1191/M1192



5 µl/lane, 8 cm gel
1×TAE, 7 V/cm, 45min

Advised Electrophoresis Conditions

Specifications

- ▶ Mix Concentration: 144ng/µl
- ▶ Recommended Loading: 2-5µl/Lane
- ▶ Content (bp): 100, 200, 300, 400, 500, 700, 1000, 1500, 2000, 3000, 4000, 5000, 6000, 8000, 10000
- ▶ Concentration: Typical Bands 100ng/5µl & Other Bands 40ng/5µl
- ▶ Storage: Stable for 3 months at RT, for long-time storage, please store at -20°C.

CAT: M1071 or M1191

**Must purchase any 10 units of 50ug GB-Ruler™ 100 bp DNA Ladder Plus or GB-Ruler™ 1 kb DNA Ladder Plus (CAT: M1071, M1191) at full price for the promotion to apply. Offer valid until December 9, 2022*



GeneBio
Systems, Inc.

Email order@genebiosystems.com
with proof of purchase for **5 FREE**

DNA Ladders